2A2, 2A3, 2C2, 2C3

INSTRUCTION FOR USE (IFU)

EN 2

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INFORMATION

Instructions for use last updated: 2026-01-03

- ▶ Read this document carefully before using the product.
- ▶ To avoid injury and damage to the product, the safety instructions must be followed.
- ▶ Instruct the user on the correct and safe use of the product.
- ▶ Save this document.

The product "tube adapter 2A2/2A3/2C2/2C3)" is hereinafter simply called the product/tube adapter.

These instructions for use contain important information regarding the application, adaptation and handling of the product.

Use the product only in accordance with the information contained in the accompanying documents.

2 Product Description

2.1 Function

Tube adapters are used as structural elements of modular prostheses of the lower extremities. They are connected by modular connecting elements.

The combination of modules allows you to control angular movements at the sagittal and frontal levels, as well as transfer loads between modules.

These instructions for use apply to the following types of tube adapters:

1. 2A2, 2A3;
2. 2C2, 2C3.

2.2 Combination capabilities

This prosthetic component is compatible with Parashar Industries' modular component system.

Compatibility with components from other manufacturers that have compatible modular connectors has not been tested.

3 Intended Use

3.1 Indications for use

The product is intended exclusively for exoprosthetics of the lower extremities.

3.2 Terms of Use

The product is designed for everyday use and should not be used for unusual activities (e.g. extreme sports: rock climbing, parachuting, paragliding, etc.).

Permissible environmental conditions are described in the technical data.

Our components work optimally in combination with the corresponding components according to body weight and mobility level according to the K-level classification, which have suitable modular connectors.

The product is recommended for mobility levels 1 (indoor walking), 2 (restricted outdoor walking), 3 (outdoor free walking) and 4 (outdoor free walking with particularly high requirements).

The allowed maximum body weight is 125 kg.

3.3 Indications

- For patients with amputation of the lower leg, thigh, dismemberment in the hip joint;
- For unilateral or bilateral amputation.

3.4 Contraindications

- All conditions that contradict or exceed the specifications specified in the "Safety" and "Indications for Use" sections.

3.5 Qualifications

A patient can only be tried on with the product by a certified specialist in orthotics and prosthetics.

4 Safety

4.1 Explanation of Warning Symbols

Symbol	Meaning
WARNINGS	Warning of possible serious risks of accident or injury.
CAUTION	Warning of possible risks of accident or injury.
WARNING	Warning of possible technical damage.

4.2 Safety Instruction Structure

WARNINGS

The title describes the source and/or type of hazard. The introduction describes the consequences of not following safety instructions.

4.3 General Safety Instructions

WARNING

Wrong combination of components

- Risk of injury due to improper combination of components and selection of individual parts.
- Use the product only according to the application.
- Follow the instructions for use.
- Combine components only within acceptable combinations.

WARNINGS

- Risk of injury and risk of damage to the product.
- Observe the conditions of use of the product and do not subject it to excessive stress.
- Follow the combination/combo exclusion instructions specified in the instructions for use of the product.
- Do not use the product in unacceptable environmental conditions.

- ▶ Check the product for damage if it has been subjected to unacceptable conditions.
- ▶ Do not use the product if it is damaged or in questionable condition. Take appropriate measures (cleaning, repair, replacement, inspection by the manufacturer or specialist workshop).
- ▶ Observe the maximum service life of the product.
- ▶ Handle the product with care to avoid mechanical damage.
- ▶ If damage to the product is suspected, check its performance.
- ▶ Do not use the product with limited performance.

Signs of change or loss of function

Changes in function can manifest themselves, for example, in the form of a change in gait, the location of the prosthesis components relative to each other, as well as in the form of noise.

WARNINGS

Fitting in the presence of children

Trauma due to swallowing small parts.

- ▶ Never leave children unattended during the fitting.

WARNING

Self-intervention in the product

Falling due to failure of loaded components or product failure.

- ▶ Any manipulation of the product, except for those described in these instructions, is prohibited;
- ▶ Only the Parashar Industries service center has the right to disassemble the product.

4.4 Patient Information

WARNING

Penetration of dirt and moisture into the product

Falling due to unexpected product behavior due to malfunction.

- ▶ Do not allow solid particles or liquids to enter the product.

WARNING**Overload due to unusual activities**

Fall due to failure of loaded components.

- ▶ The product is designed for everyday use and should not be used for extreme activities.
- ▶ Careful handling of the product and its components not only prolongs the service life, but also ensures your personal safety!
- ▶ If the product has been subjected to extreme loads (e.g. falls), it must be checked immediately for damage by a prosthetist. If necessary, the product will be sent to the Parashar Industries service center.

WARNING**Improper product care**

Damage to the product due to the use of improper cleaning products and/or disinfectants.

- ▶ Clean the product only with a damp cloth and mild soap.
- ▶ Disinfect the product with isopropyl alcohol as needed.
- ▶ Do not use disinfectants based on hydrogen peroxide and silver salts.
- ▶ Use only fresh water with a temperature not exceeding 65 °C for cleaning.

5 Package Contents and Accessories**5.1 Package Contents**

1. 1 pc. tube adapter 2A2/2A3/2C2/2C3;
2. 4 pcs. set screws M-8x16 to 2A2/2A3 and M-8x12 to 2C2/2C3;
3. 1 pc. Instructions for use (for qualified personnel).

5.2 Accessories (supplied separately).**5.3 Spare parts (supplied separately).**

6 Preparing the Product for Use

6.1 Tools and Materials Required

1. Torque wrench up to 30 Nm;
2. Pipe cutter;
3. Chamfer d=30;
4. Adhesive Loctite 243.

6.2 Mounting the Carrier Module

6.2.1 Preparing for Installation

Shorten the tube to the desired length.

Deburring the outside and inside of the edge with a deburring device.

WARNING

Improper assembly or installation

- Risk of injury due to defective prosthesis components.
- Pay attention to the installation instructions.

WARNING

Incorrect installation of threaded connections

- Risk of injury due to rupture or loose threaded connections.
- The threads should be cleaned each time before installation.
- Observe the moments when the settings are tightened.
- Pay attention to the instructions for screw length and preservation of threaded connections.

WARNING

Improper processing of the tube adapter

- Fall due to damage to the tube adapter.
- Do not clamp the tube adapter in a vice during processing.
- The tube should only be shortened with a pipe cutter or a special device.

6.2.2 Installation in a modular prosthesis

WARNING

Incorrect installation in a modular prosthesis
Risk of injury due to destruction of load-bearing parts
During installation, the carrier module should be pushed all the way into the prosthetic component designed for this purpose

1. Place the socket for the alignment pyramid of the carrier module in the distal part of the prosthesis
2. Align the screw module clamping groove:
Screw module: in the front plane
Screw module, movable: in the anterior or medial plane
3. Connect the components to the carrier module according to the instructions.

Connecting the Alignment Pyramid to the Socket for the Alignment Pyramid

Required materials:
torque wrench), Loctite 243 glue

1) Fitting:

Screw in the threaded screws.

1. Tighten the threaded screws with torque wrench:(10Nm) for 2C2,2C3; (15 Nm) for 2A2,2A3

2) Final installation:

Apply Loctite 243 glue to the threaded screws.

Screw in the threaded screws.

Pre-tighten the threaded screws with a torque wrench (10Nm), then tighten tightly (12Nm) for 2C2,2C3;

Pre-tighten the threaded screws with a torque wrench (15Nm), then tighten tightly (18Nm) for 2A2,2A3.

Replacement and dismantling

During replacement or disassembly, the corrected position of the prosthesis component can be maintained. To do this, unscrew both threaded screws screwed deeper than the others and located next to each other.

7 Maintenance

Annual safety inspections are recommended to ensure patient safety and product performance.

8 Cleaning

1. Clean the product with a damp cloth and soap.
2. Rinse with clean fresh water.
3. Dry with a soft cloth.
4. **NOTE!** Do not store the product wet. Allow it to air dry.

9 Disposal

In some countries, it is forbidden to dispose of these products with household waste. Improper disposal can be harmful to health and the environment. Follow the instructions of the national authority for return and collection.

10 Legal Information

All legal terms are subject to the laws of the country of use.

10.1 Liability

The manufacturer is only responsible if the product is used in accordance with these instructions.

The manufacturer is not responsible for damage caused by non-compliance with the instructions, improper use or unauthorized intervention.

10.2 Trademarks

All product names are trademarks owned by Parashar Industries.

10.3 CE Compliance

Parashar Industries declares that the product complies with current European requirements for medical devices.

10.4 Compliance with Technical Regulations

Parashar Industries guarantees compliance with the Technical Regulations on Medical Devices No. 753 dated October 2, 2013, provided that the consumer complies with the conditions of transportation, storage, installation, and operation.

10.5 Warranty

The manufacturer provides a warranty for this product from the date of purchase.

The warranty covers failures caused by confirmed defects in material, manufacturing defects or design defects that have been notified to the manufacturer during the warranty period. Detailed information on the terms of the warranty can be obtained from the relevant representative of the manufacturer or from the company supplying the products.

11 Technical Data

Environmental conditions

Permissible conditions for using the product:

Temperature range of use from -10°C to +60°C

Storage and transportation conditions: from -10°C to +60°C

Fresh water, soapy water, urine (needs to be washed)

Swimming pool (needs flushing)

Sea and chlorinated water (requires flushing)

Unacceptable conditions for using the product:

Mechanical vibration or mechanical shock

Ingress of dust, sand, hygroscopic parts (e.g. talc)

Acidic and alkaline solutions

Solvent-containing cleaners

Product characteristics

Parameter	2A2	2A3	2C2	2C3
Weight [g]	213	344	222	354
Proximal connection	Pyramid Receiver/Pipe d=30			
Tightening torque of the screws of the proximal connection [N·m]	18		12	
Distal junction	Pipe d=30/ Pyramid receiver			
Max. Body weight [kg]	125			
System Height [mm]	235	435	233	433
Structure height [mm]				
Tube material	Aluminum alloy			
Adapter Material	Aluminum alloy		Steel	
Service life [years]	3			

12 Appendices

12.1 Symbols Used

 LOT	Batch number		Keep out of direct sunlight
	Manufacturer		Protect from moisture
	Date of manufacture		Medical device
	See instructions for use		Conformity mark without an authority code
	Multiple uses, but for one patient		Keep out of direct sunlight
	Declaration of conformity to European directives		

Parashar Industries LLC

12 Makarivska Street, Petropavlivska
Borshchahivka village, Bucha District,
Kyiv Region, 08129, Ukraine

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