EAGLE-X


INSTRUCTION FOR USE (IFU)


EN 5

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






INFORMATION

Instructions for use last updated: 2026-01-03

- Read this document carefully before using the product and follow the safety warnings.
- Teach the user how to use the product safely.
- Contact the manufacturer if you have any questions about the product or have any problems.
- Report every serious incident related to the product to the manufacturer and the relevant authority in your country. This is especially important when health is deteriorating.
- Save this document for your records.

The product " Knee joint EAGLE-X " is hereinafter simply called the product/knee joint.

These instructions for use contain important information regarding the application, adaptation and handling of the product.

Use the product only in accordance with the information contained in the accompanying documents.

2 Product Description

2.1 Function

The EAGLE-X knee joint is a polycentric modular knee assembly with a lock, rotary hydraulic system, geometric lock and SES (Smooth Ergonomic Stance) cushioning function, which is used for prosthetics after amputation above the knee.

A hydraulic system controls the transfer phase. The level of damping of flexion, extension and support is adjusted individually. This ensures a comfortable ride over a wide range of speeds. The softness of the shock absorption can be changed by adjusting the SES shock absorber block. The EAGLE-X knee assembly implements the concept of a geometric lock to increase the level of strut resistance in the support phase.

The EAGLE-X knee joint is designed to meet the increased requirements for prosthetics, as it has a greater range of adjustments and, accordingly, is able to provide more dynamic and physiological walking with less energy consumption for the user.

2.2 Combination capabilities

This prosthetic component is compatible with Parashar Industries' modular component system. Compatibility with components from other manufacturers that have compatible modular connectors has not been tested.

3 Intended Use

3.1 Indications for use

The product is intended exclusively for exoprosthetics of the lower extremities.

3.2 Terms of Use

The product is designed for everyday use and should not be used for unusual activities (e.g. extreme sports: rock climbing, parachuting, paragliding, etc.).

ATTENTION!

Product overload

- Use the product only in accordance with the conditions of use specified by the manufacturer

Permissible conditions for using the product:

Temperature range of use from -10°C to +60°C

Storage and transportation conditions: from -10°C to +60°C

Fresh water, soapy water, urine (needs to be washed)

Swimming pool (needs flushing)

Sea and chlorinated water (requires flushing)

Unacceptable conditions for using the product:

Mechanical vibration or mechanical shock

Ingress of dust, sand, hygroscopic parts (e.g. talc)

Acidic and alkaline solutions

Solvent-containing cleaners

Our components work optimally in combination with the corresponding components according to body weight and mobility level according to the K-level classification, which have suitable modular connectors.

The product is recommended for mobility level 2 (limited outdoor walking), 3 (free outdoor walking).

The allowed maximum body weight is 125 kg.

3.3 Indications

- For patients with hip amputation, dismemberment in the hip joint;
- For unilateral or bilateral amputation.

3.4 Contraindications

- All conditions that contradict or exceed the specifications specified in the "Safety" and "Indications for Use" sections.

3.5 Qualifications

A patient can only be tried on with the product by a certified specialist in orthotics and prosthetics.

4 Safety**4.1 Explanation of Warning Symbols**

Symbol	Meaning
WARNINGS	Warning of possible serious risks of accident or injury.
WARNING	Warning of possible risks of accident or injury.
REMARKS	Warning of possible technical damage.

4.2 Safety Instruction Structure**WARNINGS**

The title describes the source and/or type of hazard. The introduction describes the consequences of not following safety instructions.

4.3 General Safety Instructions**WARNINGS****Failure to comply with safety warnings**

Injury/damage to the product due to the use of the product in certain situations.

- Follow the safety warnings and precautions specified in this document.

WARNINGS**Fitting in the presence of children**

Trauma due to swallowing small parts.

- Never leave children unattended during the fitting.

WARNING

Self-intervention in the product

Falling due to failure of loaded components or product failure.

- Any manipulation of the product, except for those described in these instructions, is prohibited;
- Only the Parashar Industries service center has the right to disassemble the product.

WARNING

Reuse of the product by another patient

Falls due to loss of function and damage to the product.

- The product is only allowed to be used by one patient.

WARNING

Exceeding the service life

Falls due to alteration or loss of function, as well as damage to the product.

- Care should be taken not to exceed the service life specified in this section.

WARNING

Clamping in the area of the knee mechanism

Pinching of limbs (e.g. fingers) and skin due to uncontrolled movement of the product.

- Never handle the product in the hinge area during use.
- Installation and adjustment should only be carried out in conditions of increased care.

WARNING

Mechanical damage to the product

Risk of injury due to alteration or loss of function.

- Handle the product with care.
- Check the damaged product for operability and the possibility of further use.
- Do not use the product when changing or losing functions (see Fig. "Signs of change or loss of functions during operation" in this section).
- If necessary, take appropriate measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of change or loss of function during operation

Changes in the operation of the product can be manifested, for example, in a change in gait (transfer phase, support phase), incomplete extension, tightness of gait and the appearance of noises.

WARNING

Use of incompatible prosthetic components

Changes in the operation of the product can be manifested, for example, in a change in gait (transfer phase, support phase), incomplete extension, tightness of gait and the appearance of noises.

4.4 Patient Information

WARNING

Penetration of dirt and moisture into the product

Falling due to unexpected product behavior due to malfunction.

- Do not allow solid particles or liquids to enter the product.

WARNING

Overload due to unusual activities

Fall due to failure of loaded components.

- The product is designed for everyday use and should not be used for extreme activities.
- Careful handling of the product and its components not only prolongs the service life, but also ensures your personal safety!
- If the product has been subjected to extreme loads (e.g. falls), it must be checked immediately for damage by a prosthetist. If necessary, the product will be sent to the Parashar Industries service center.

REMARKS

Improper product care

Damage to the product due to the use of improper cleaning products and/or disinfectants.

- Clean the product only with a damp cloth and mild soap.
- Disinfect the product with isopropyl alcohol as needed.
- Do not use disinfectants based on hydrogen peroxide and silver salts.
- Use only fresh water with a temperature not exceeding 65 °C for cleaning.

5 Package Contents and Accessories

5.1 Package Contents

1. 1 pc. knee joint Eagle-X;
2. 1 pc. Instructions for use (for qualified personnel).

5.2 Accessories (supplied separately).

5.3 Spare parts (supplied separately).

6 Preparing the Product for Use

6.1 Required tools

1. Torque wrench up to 30 Nm;
2. Adhesive Loctite 243.

6.2 Installation

6.2.1 Prosthesis Manufacturing Guidelines

WARNING

Improper assembly, installation or adjustment

- Injuries due to improper installation, adjustment or damage to the components of the prosthesis.
- Pay attention to the installation, installation and adjustment instructions.

WARNING

First use by the patient

- Falls due to insufficient patient experience or improper assembly or adjustment of the prosthesis.
- For the safety of the patient during the first attempts to stand and walk, assistive devices (e.g. walking bars, handrails or rollator) must be used.

6.2.2 Performing basic prosthesis assembly

- To securely place the patient in a standing position, assemble the prosthesis in accordance with the assembly figures and instructions for use of all components of the prosthesis.
- When activating and deactivating the lock, the knee joint must be pressed in the direction of the extension limiter.
- Perform the general layout of the prosthesis with the knee joint blocked.

- Information: The back of the knee joint can be used as a flexion limiter for the stump sleeve.
- After installing the sleeve, make sure that no metal parts are pressing on the back of the knee joint.

WARNING

Improper installation of threaded connections

Risk of injury due to breakage or loosening of threaded connections.

- Clean the threads each time before installation.
- Observe the set tightening torques.
- Pay attention to the instructions for fixing threaded connections.

WARNING

Improper processing of the carrier module

Fall due to damage to the carrier module.

- During processing, do not clamp the load-bearing module in a vice.
- Shorten the load-bearing module only using a pipe cutter or a special cutting device.
- Deburring the outside and inside of the cut edge using a deburring tool.

WARNING

Improper installation of the carrier module

Risk of injury due to the destruction of load-bearing parts.

- Clean the contact surfaces of the carrier module and the tube mount with a degreasing agent.

6.2.2 Basic assembly

6.2.3 Basic assembly

Fig. 1 illustrates the maximum bend of the knee joint.

In the construction of a prosthesis of the lower limb using the EAGLE-X knee node (Fig. 2), it is not recommended to use a cosmetic coating (polyurethane foam).

6.2.4 EAGLE link mechanism

During the support phase, at the very beginning of the knee load, immediately after the heel comes into contact with the floor, the shock absorption function affects the dynamics of the node (Fig. 3a). This is due to the fact that the instantaneous center of rotation (the point of intersection of the lines extending the links is marked in the figures with dotted lines) moves even more beyond the load line and thus random bending is blocked by the geometric lock even more, which makes the support phase safer.

The support reaction force that occurs when the heel comes into contact with the floor moves the SES shock absorption unit forward, causing the center of instantaneous rotation to be moved up and back further away from the load line (see 3b).

When the heel is detached, the SES shock absorption unit, together with the instantaneous center of rotation, returns back to its original position (Fig. 4a), which allows flexion of the knee joint before the start of the transfer phase (Fig. 4b).

Geometric lock

Due to its design, the EAGLE-X knee assembly implements the concept of a geometric lock to increase the level of strut resistance in the support phase.

When the heel comes into contact, the lower part of the knee joint bends slightly (Fig. 4d), which minimizes the shock load on the stump. As a result, the instantaneous center of rotation of the knee joint moves up and back from the load line (Fig. 4c), which further prevents the likelihood of knee flexion at the beginning of the support phase.

Inside the support phase, the SES cushioning unit returns to its original position, allowing the knee to be flexed before the start of the transfer phase (Fig. 4c).

SES (Smooth Ergonomic Stance) damping unit function

Starting from the moment the heel comes into contact with the floor and moves to the middle position of the support phase, the SES shock absorption unit is activated, which absorbs the impact of the foot while walking, maintaining a smooth forward movement of the body and has a natural effect, minimizing energy costs. The softness of the shock absorption can be changed by adjusting the SES shock absorber block.

Hydraulic cylinder

The compact hydraulic cylinder provides a facilitated transfer phase. The hydraulic cylinder is designed to ensure smooth walking in the transfer phase and shows good adaptability to different walking speeds.

6.2.5 Bench adjustment

To take full advantage of the EAGLE-X knee joint, it is imperative to properly configure it for the user. We recommend the following bench adjustment scheme:

1. In the sagittal plane (side view).

The load line (a sheer line lowered from the midpoint of the inner part of the receiver stump) should pass through the center of the upper axis of the anterior vertical link (Fig. 5a), reaching the middle of the imaginary line between the end of the heel and the ball of the foot (Fig. 5b).

2. In the frontal plane (rear view).

The load line (a sheer line running from the midpoint of the inner part of the stump sleeve) should pass through the center of the knee and the center of the foot (Fig. 5b).

6.2.6 Dynamic fitting

First you need to adjust the level of strut resistance of the knee assembly in the support phase.

The strut resistance of a fully extended knee joint at the time of heel impact significantly depends on the inclination of the axis of the knee joint.

In the case of the EAGLE-X knee joint, the instantaneous center of rotation significantly changes position immediately after the prosthesis is loaded, which makes the knee joint safe. A series of test walks is carried out, including adjustment of the product, looking at the features that may be detected.

For quick and efficient adjustment of the EAGLE-X knee joint, a number of step-by-step adjustments are recommended.

STEP 1: Adjusting the strut stability by turning the knee node in the sagittal plane.

The first to adjust the strut stability is by turning the knee node in the sagittal plane.

During a test walk, if you notice such negative manifestations of the knee node as:

1. too much depreciation, accompanied by mechanical noise;
2. difficulty bending the knee at the moment of lifting the toe off the floor.

To eliminate negative manifestations, it is necessary to reduce the strut resistance by turning the knee joint clockwise (Fig. 6a). In this case, the instantaneous center of rotation will shift forward.

If during a test walk you notice such negative manifestations of the knee node as:

1. the user feels instability at the moment of impact of the heel with the floor;
2. There is no appropriate cushioning when walking on an inclined surface.

To eliminate negative manifestations, it is necessary to increase the strut resistance by turning the knee joint counterclockwise (Fig. 6b). In this case, the instantaneous center of rotation will shift backwards.

STEP 2: Adjusting the stiffness of the SES suspension unit.

During a test walk, the patient may notice that the feeling of posture is not rigid enough (reliable). The feeling of instability is usually caused by bouncing movements of the knee knot. If the forward tilt of the knee joint is already sufficient, and there is no risk of accidental bending of the knot, then it is the increase in the hardness of the SES shock absorber block (increase in the initial compression of the shock absorber) that can help eliminate discomfort (Fig. 7).

To do this, follow the instructions for setting up the knee joint:

1. Turn the M8 screw clockwise (located behind the knee assembly) 1/4 turn (90 degrees) using a 4 mm wrench.
2. If the result is not enough, repeat this procedure.

STEP 3. Adjusting the hydraulic unit.

Adjustment of the hydraulic resistance of the EAGLE-X knee joint is possible separately for flexion and extension.

Extension adjustment controls the rate of extension of the knee module. After maximum flexion of the node in the transfer phase, the knee module is quickly straightened and gently dampened before full straightening, avoiding terminal impact (Fig. 8).

Flexion adjustment controls the angle of flexion of the prosthesis when walking (Fig. 9).

ATTENTION! Clockwise rotation of the adjustment screws increases the flexion and extension resistance, and counterclockwise reduces the flexion and extension resistance accordingly (Fig.8/Fig.9)

STEP 4. Fixation of the EAGLE-X knee joint with a locking mechanism

The EAGLE-X knee module is fixed by switching the lock lever to the Lock position (Fig. 10). The lock is defixed in the Unlock position (Fig. 11).

7 Maintenance

Annual safety inspections are recommended to ensure patient safety and product performance.

8 Cleaning

1. Clean the product with a damp cloth and soap.
2. Rinse with clean fresh water.
3. Dry with a soft cloth.
4. **NOTE!** Do not store the product wet. Allow it to air dry.

9 Disposal

In some countries, it is forbidden to dispose of these products with household waste. Improper disposal can be harmful to health and the environment. Follow the instructions of the national authority for return and collection.

10 Legal Information

All legal terms are subject to the laws of the country of use.

10.1 Liability

The manufacturer is only responsible if the product is used in accordance with these instructions.

The manufacturer is not responsible for damage caused by non-compliance with the instructions, improper use, combination with components from other manufacturers, or unauthorized intervention.

10.2 Trademarks

All product names are trademarks owned by Parashar Industries.

10.3 CE Compliance

Parashar Industries declares that the product complies with current European requirements for medical devices.

10.4 Compliance with the Technical Regulations of Ukraine

The manufacturer Parashar Industries confirms that the product meets the requirements of the Technical Regulation on Medical Devices, approved by the Resolution of the Cabinet of Ministers of Ukraine No. 753 dated October 2, 2013, provided that the user complies with the established requirements for transportation, storage, installation and operation of the product.

10.5 Warranty

The manufacturer provides a warranty for this product from the date of purchase.

The warranty covers failures caused by confirmed defects in material, manufacturing defects or design defects that have been notified to the manufacturer during the warranty period.

Detailed information on the terms of the warranty can be obtained from the relevant representative of the manufacturer or from the company supplying the products.

11 Technical Data

Environmental conditions

Product characteristics

Parameter	Knee joint EAGLE-X
Weight [g]	990
Proximal connection	Aligning pyramid
Distal junction	Aligning pyramid
Max. Body weight [kg]	125
Mobility level	K2, K3
System Height [mm]	145,5
System height in the proximal part to the starting point of assembly [mm]	3
System height in the distal part to the starting point of assembly [mm]	142,5
Maximum bending angle	135°
Material	Aluminum alloy
Service life [years]	3

12 Appendices

12.1 Symbols Used

SN Batch number	 Keep out of direct sunlight
 Manufacturer	 Protect from moisture
 Date of manufacture	MD Medical device
 See instructions for use	 Conformity mark without an authority code
 Multiple uses, but for one patient	 Keep out of direct sunlight
CE Declaration of conformity to European directives	

Parashar Industries LLC
12 Makarivska Street, Petropavlivska
Borshchahivka village, Bucha District,
Kyiv Region, 08129, Ukraine
tel.: +38 (098) 227 2659
info@parasind.com
www.parasind.com

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