PARASHAR

MAYA PK-72






INSTRUCTION FOR USE (IFU)

EN 4

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
UKR 12

4

INFORMATION**Instructions for use last updated: 2026-01-03**

- Read this document carefully before using the product and follow the safety warnings.
- Teach the user how to use the product safely.
- Contact the manufacturer if you have any questions about the product or have any problems.
- Report every serious incident related to the product to the manufacturer and the relevant authority in your country. This is especially important when health is deteriorating.
- Save this document for your records.

The product "knee joint MAYA PK-72" is hereinafter simply called the product/knee joint.

These instructions for use contain important information regarding the application, adaptation and handling of the product.

Use the product only in accordance with the information contained in the accompanying documents.

2 Product Description**2.1 Function**

The MAYA PK-72 knee joint is a monocentric knee joint with a hydraulic cylinder and a mechanical lock, which is used for prosthetics after amputation above the knee.

The hydraulic cylinder of the knee module increases the safety of use and simulates the natural seating mechanism. The maximum extension angle of the knee module is 130° Fig.1 s.

The use of aluminum and stainless steel in the production made it possible to reduce the weight of the knee joint and made it less whimsical to maintain.

2.2 Combination capabilities

This prosthetic component is compatible with Parashar Industries' modular component system. Compatibility with components from other manufacturers that have compatible modular connectors has not been tested.

3 Intended Use**3.1 Indications for use**

The product is intended exclusively for exoprosthetics of the lower extremities.

3.2 Terms of Use

The product is designed for everyday use and should not be used for unusual activities (e.g. extreme sports: rock climbing, parachuting, paragliding, etc.).

ATTENTION!**Product overload**

- Use the product only in accordance with the conditions of use specified by the manufacturer

Permissible conditions for using the product:

Temperature range of use from -10°C to +60°C

Storage and transportation conditions: from -10°C to +60°C

Fresh water, soapy water, urine (needs to be washed)

Swimming pool (needs flushing)

Sea and chlorinated water (requires flushing)

Unacceptable conditions for using the product:

Mechanical vibration or mechanical shock

Ingress of dust, sand, hygroscopic parts (e.g. talc)

Acidic and alkaline solutions

Solvent-containing cleaners

Our components work optimally in combination with the corresponding components according to body weight and mobility level according to the K-level classification, which have suitable modular connectors.

The product is recommended for mobility levels 1 (indoor walking), 2 (restricted outdoor walking), 3 (outdoor free walking).

The allowed maximum body weight is 125 kg.

3.3 Indications

- For patients with hip amputation, dismemberment in the hip joint;
- For unilateral or bilateral amputation.

3.4 Contraindications

- All conditions that contradict or exceed the specifications specified in the "Safety" and "Indications for Use" sections.

3.5 Qualifications

A patient can only be tried on with the product by a certified specialist in orthotics and prosthetics.

4 Safety

4.1 Explanation of Warning Symbols

Symbol	Meaning
WARNINGS	Warning of possible serious risks of accident or injury.
WARNING	Warning of possible risks of accident or injury.
REMARKS	Warning of possible technical damage.

4.2 Safety Instruction Structure

WARNINGS

The title describes the source and/or type of hazard. The introduction describes the consequences of not following safety instructions.

4.3 General Safety Instructions

WARNINGS

Failure to comply with safety warnings

Injury/damage to the product due to the use of the product in certain situations.

- Follow the safety warnings and precautions specified in this document.

WARNINGS

Fitting in the presence of children

Trauma due to swallowing small parts.

- Never leave children unattended during the fitting.

WARNING

Self-intervention in the product

Falling due to failure of loaded components or product failure.

- Any manipulation of the product, except for those described in these instructions, is prohibited;
- Only the Parashar Industries service center has the right to disassemble the product.

WARNING

Reuse of the product by another patient

Falls due to loss of function and damage to the product.

- The product is only allowed to be used by one patient.

WARNING

Exceeding the service life

Falls due to alteration or loss of function, as well as damage to the product.

- Care should be taken not to exceed the service life specified in this section.

WARNING

Clamping in the area of the knee mechanism

Pinching of limbs (e.g. fingers) and skin due to uncontrolled movement of the product.

- Never handle the product in the hinge area during use.
- Installation and adjustment should only be carried out in conditions of increased care.

WARNING

Mechanical damage to the product

Risk of injury due to alteration or loss of function.

- Handle the product with care.
- Check the damaged product for operability and the possibility of further use.
- Do not use the product when changing or losing functions (see Fig. "Signs of change or loss of functions during operation" in this section).
- If necessary, take appropriate measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of change or loss of function during operation

Changes in the operation of the product can be manifested, for example, in a change in gait (transfer phase, support phase), incomplete extension, tightness of gait and the appearance of noises.

WARNING

Use of incompatible prosthetic components

Falls due to unexpected product behavior or failure of loaded components.

- Only use the product in combination with components approved by Parashar Industries (see instructions and catalogs).

4.4 Patient Information

WARNING

Penetration of dirt and moisture into the product

Falling due to unexpected product behavior due to malfunction.

- Do not allow solid particles or liquids to enter the product.

WARNING

Overload due to unusual activities

Fall due to failure of loaded components.

- The product is designed for everyday use and should not be used for extreme activities.
- Careful handling of the product and its components not only prolongs the service life, but also ensures your personal safety!
- If the product has been subjected to extreme loads (e.g. falls), it must be checked immediately for damage by a prosthetist. If necessary, the product will be sent to the Parashar Industries service center.

REMARKS

Improper product care

Damage to the product due to the use of improper cleaning products and/or disinfectants.

- Clean the product only with a damp cloth and mild soap.
- Disinfect the product with isopropyl alcohol as needed.
- Do not use disinfectants based on hydrogen peroxide and silver salts.
- Use only fresh water with a temperature not exceeding 65 °C for cleaning.

5 Package Contents and Accessories

5.1 Package Contents

1. 1 pc. knee joint MAYA PK-72;
2. 1 pc sweeping system;
3. 1 pc. traction cable;
4. 1 pc. Instructions for use (for qualified personnel).

5.2 Accessories (supplied separately).

5.3 Spare parts (supplied separately).

6 Preparing the Product for Use

6.1 Required tools

1. Torque wrench up to 30 Nm;
2. Adhesive Loctite 243.

6.2 Installation

6.2.1 Prosthesis Manufacturing Guidelines

WARNING

Improper assembly, installation or adjustment

- Injuries due to improper installation, adjustment or damage to the components of the prosthesis.
- Pay attention to the installation, installation and adjustment instructions.

WARNING

First use by the patient

- Falls due to insufficient patient experience or improper assembly or adjustment of the prosthesis.
- For the safety of the patient during the first attempts to stand and walk, assistive devices (e.g. walking bars, handrails or rollator) must be used.

6.2.2 Performing basic prosthesis assembly

- To securely place the patient in a standing position, assemble the prosthesis in accordance with the assembly figures and instructions for use of all components of the prosthesis.
- When activating and deactivating the lock, the knee joint must be pressed in the direction of the extension limiter.

- Perform the general layout of the prosthesis with the knee joint blocked.
- Information: The back of the knee joint can be used as a flexion limiter for the stump sleeve.
- After installing the sleeve, make sure that no metal parts are pressing on the back of the knee joint.

WARNING

Improper installation of threaded connections

Risk of injury due to breakage or loosening of threaded connections.

- Clean the threads each time before installation.
- Observe the set tightening torques.
- Pay attention to the instructions for fixing threaded connections.

WARNING

Improper processing of the carrier module

Fall due to damage to the carrier module.

- During processing, do not clamp the load-bearing module in a vice.
- Shorten the load-bearing module only using a pipe cutter or a special cutting device.
- Deburring the outside and inside of the cut edge using a deburring tool.

WARNING

Improper installation of the carrier module

Risk of injury due to the destruction of load-bearing parts.

- Clean the contact surfaces of the carrier module and the tube mount with a degreasing agent.

6.2.2 Guidelines for the use of cosmetics

In the construction of a prosthesis of the lower limb using the MAYA PK-72 knee node, it is not recommended to use a cosmetic coating (polyurethane foam).

6.2.3 Basic assembly

The spatial arrangement of the stump sleeve of the prosthesis and modular components affects the static and dynamic characteristics of the prosthesis. The position of the axis of the knee node affects its characteristics.

Therefore, please note that in the extended state, the load line is in front of the axis of the knee joint, which ensures sufficient stability and the possibility of unblocking by applying a slight force.

When determining the position of the stump sleeve, the position of the stump of the lower limb should be taken into account.

The assembly of the prosthesis and its fit, as a rule, are performed in the following sequence:

1. Basic assembly;
2. Bench adjustment;
3. Dynamic fitting.

In the case of primary prosthetics, the above stages are performed with a trial (training) prosthesis, for early verticalization and ensuring the best construction to ensure a natural gait. After fitting the stump sleeve and ensuring the necessary adjustment of the prosthesis components, the main prosthesis is made for the user. In this case, all manufacturing steps are repeated.

6.2.4 Knee assembly installation

The proximal connection (with the adapter) is made using the alignment screws. Distal - using a clamp for a pipe $d=30$ and an M6 screw.

The structure of the MAYA PK-72 knee joint is shown in Fig. 2:

1. proximal connection;
2. lock control lever;
3. distal connection;
4. Clamping screw.

Clamping of the screw on the elbow assembly is carried out exclusively with a torque wrench and must be taken into account the manufacturer's marking, which is a caveat, regarding the recommended clamping force of 10 Nm.

6.2.5. Configuring the "permanent lock" mode

In the "Permanent Lock" mode, the knee joint has a mechanical flexion lock built into the upper part of the product, which automatically holds the product in an extended position (Fig. 1 a), while its closure is noticeable for the patient.

The flexion is unlocked by a traction cable attached to the lever (Fig. 2 p. 4) and is carried out when it is pulled.

The user has the option of unlocking the knee joint before sitting down, or already in a sitting position. To unlock the knee node, the patient must almost completely unload the prosthesis by transferring weight to the opposite limb, followed by flexion of the stump. When unlocked, the user can perform a bend. In a standing position, the user automatically locks the module by applying an extension force with the prosthesis fully extended.

6.2.6. Squat mode

A feature of the MAYA PK-72 knee joint is the presence of a hydraulic cylinder that acts as an auxiliary mechanism for a natural squat. Depending on the anthropometric data of the user, it is possible to individually adjust the resistance of the hydraulic cylinder (squat speed).

To adjust the resistance of the hydraulic cylinder, use the regulator (Fig. 3a p.1).

Ask the patient to open the knee lock and sit down.

If the resistance is too small, turn the regulator clockwise;

If the resistance is too high, turn the regulator counterclockwise.

6.2.7. Installation of the lock lever system

When supplying the MAYA PK-72 knee joint, it has a traction cable that must be cut and attached to the lock control lever, depending on the side of use. Shorten the cable exactly enough so that it is slightly taut, but does not unlock the knee lock. The lock should be opened only after pulling the cable at a distance of about 5 mm.

The traction cable must be located on the lateral side of the stump sleeve.

6.2.8 Attaching the traction cable to the knee assembly

Pass the pull cable through the hole in the lever (Fig. 4a).

Install the bushing on the cable (Fig. 4a).

Attach the pull cable by warming it up and pressing it with its end on a flat surface (Fig. 4b).

6.2.9 Assembling the Latch System

Steps for assembling the lock lever system:

1. Install the lock part (Fig. 5, item 4) into the stump sleeve when laminating;
2. Determine the length of the cable (Fig. 5, item 1);
3. Cut the cable to a suitable length;
4. Heat the cable and pass it around the screw (Fig. 5 point 3), installing it in the lock lever (Fig. 5 point 2);
5. Screw the screw (Fig. 5 point 2) into the base of the lock lever (Fig. 5 point 3), fixing the cable (Fig. 5 point 1);
6. Install (Fig. 5, item 2) on the stump receiver using a screw (Fig. 5, item 5) through the sleeve (Fig. 5, item 6), as indicated in the diagram.

7 Maintenance

Annual safety inspections are recommended to ensure patient safety and product performance.

8 Cleaning

1. Clean the product with a damp cloth and soap.
2. Rinse with clean fresh water.
3. Dry with a soft cloth.
4. **NOTE!** Do not store the product wet. Allow it to air dry.

9 Disposal

In some countries, it is forbidden to dispose of these products with household waste. Improper disposal can be harmful to health and the environment. Follow the instructions of the national authority for return and collection.

10 Legal Information

All legal terms are subject to the laws of the country of use.

10.1 Liability

The manufacturer is only responsible if the product is used in accordance with these instructions.

The manufacturer is not responsible for damage caused by non-compliance with the instructions, improper use, combination with components from other manufacturers, or unauthorized intervention.

10.2 Trademarks

All product names are trademarks owned by Parashar Industries.

10.3 CE Compliance

Parashar Industries declares that the product complies with current European requirements for medical devices.

10.4 Compliance with the Technical Regulations of Ukraine

The manufacturer Parashar Industries confirms that the product meets the requirements of the

Technical Regulation on Medical Devices, approved by the Resolution of the Cabinet of Ministers of Ukraine No. 753 dated October 2, 2013, provided that the user complies with the established requirements for transportation, storage, installation and operation of the product.

10.5 Warranty

The manufacturer provides a warranty for this product from the date of purchase.

The warranty covers failures caused by confirmed defects in material, manufacturing defects or design defects that have been notified to the manufacturer during the warranty period.

Detailed information on the terms of the warranty can be obtained from the relevant representative of the manufacturer or from the company supplying the products.

11 Technical Data

Environmental conditions

Product characteristics

Parameter	MAYA PK-72
Weight [g]	668
Proximal connection	Aligning pyramid
Distal junction	Pipe clamp d=30
Max. Body weight [kg]	125
Mobility level	K1, K2, K3
System Height [mm]	107
System height in the proximal part to the starting point of assembly [mm]	3
System height in the distal part to the starting point of assembly [mm]	104
Maximum bending angle	130°
Material	Aluminum alloy
Service life [years]	3

12 Appendices

12.1 Symbols Used

SN Batch number	 Keep out of direct sunlight
 Manufacturer	 Protect from moisture
 Date of manufacture	MD Medical device
 See instructions for use	 Conformity mark without an authority code
 Multiple uses, but for one patient	 Keep out of direct sunlight
CE Declaration of conformity to European directives	

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This image shows a full page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page, typical of notebook paper. There are no margins, text, or other markings on the page.

This image shows a single page of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]